

Tracy, CA 95304
Phone: (866) 592-0488 • FAX: (209) 834-8832
topshelforthopedics.com



The Gauntlet Pro Cinch 2.0 Ankle Stabilizer (With & W/O Hinges)

Instructions for Use

Please read these instructions completely and carefully before applying this brace.
Correct application of the brace is important for proper function.



The Gauntlet Pro Cinch 2.0 - Ankle Stabilizer


(With & W/O Hinges)

Instructions for Use


Please read these instructions completely and carefully before applying this brace.
Correct application of the brace is important for proper function.

Indications: Acute ankle sprains, chronic ankle sprains and preventive use


Application Instructions:




Step #1: Before application loosen all straps. For best performance and comfort, wear a sock between the brace and skin.




Step #2: If hinges are present, confirm placement of the footplate is centered on stay pocket of the cuff within 1/4" of the top of the sleeve on both the inside and outside near the figure 8 strap attachment. The shorter portion of the footplate is positioned towards the heel.



Step #3: Insert foot into inner sock so that the heel is fully seated. Then taking the upper posterior straps above the ankle, thread the outer strap through the slot of the other strap to secure them.




Step #4: With the ankle at a 90° angle, bring the outside (lateral) stabilizing strap forward across the top of the foot and under the heel. Secure the stabilizing strap to the opposite side of the ankle using the loops provided.



Step #5: Bring the inside (medial) stabilizing strap forward across the top of the foot and under the heel. Secure the stabilizing strap to the opposite side of the ankle using the loops provided.

Once both straps are secured, re-tension both straps at the same time while seated. Use the finger loops and pull up and away firmly on both straps, reapply straps to side of brace.



Step #6: Secure the ankle brace with the elastic cuff over the medial and lateral straps. Periodic adjustment to the tension of the medial and lateral straps may be necessary depending on individual needs.

Use & Care: Hand wash in cold water using mild soap, rinse thoroughly.

Warning: This product is to be used under the supervision of a medical professional. This device is not intended for unsupervised public use. If you experience any pain, swelling, sensation changes, or any unusual reactions while using this product, consult your medical professional immediately.

Warranty: Top Shelf Orthopedics guarantees this product to be free of defects in material and construction for a period of six months from the date of purchase.

Caution: Federal Law (U.S.A) restricts this device for sale by or on the order of a licensed healthcare professional.